

# ***SONNTEK, INC.***

## **SONNTEK P.T.C 200 PRECISION TUBING CUTTER KIT**

- *Essential tools for your Lab or Service Group*
- *Simple and Easy to Use*
- *Custom cut to the length you need, right when you need it!*

The Sonntek PTC200 is far more than a tubing cutter. A wide range of accessories enables you to grind, sharpen, sand, drill, engrave, slot and otherwise modify tubing to meet the needs of your system. The kit includes:

1. Variable speed motor (5,000-30,000 RPM, ball bearing construction)
2. Cutting collar with swivel arm and wheel alignment rod
3. 36 abrasive wheels
4. Grinding, sharpening, drilling, engraving, and other wheels
5. Soft tubing cutter
6. Centering/deburring tool
7. Diagonal cutting pliers
8. Tube straightening/bending/de-kinking tool
9. Slotted screwdriver
10. 10X Illuminated magnifier
11. Dynaseal finger-tightening universal fittings
12. Razor-blades, fittings and other accessories
13. Tubing
  - a. Stainless steel: 1/16" OD x 0.010" ID, 10 ft. (3 m)
  - b. Stainless steel: 1/16" OD x 0.020" ID, 5 ft. (1.5 m)
14. Deluxe storage case

### **PTC200 KIT**

***Call for Price (Ask for Dave)***



**SONNTEK, INC. 125 PLEASANT AVE. UPPER SADDLE RIVER, NJ 07458**  
**WWW.SONNTEK.COM PH: 201-236-9300 SALES@SONNTEK.COM**