


## Sonntek Fiber Optic Filter Test Kit

Try our new test kit tool to help troubleshoot your fiber optic spectrophotometer. Save time and money by quickly determining if you are having a flow cell problem or an issue with your detector.

Our filter test kit comes with ½” filter and everything you need to assess your HPLC/FPLC detector performance. Use without the flow cell to measure absorbance and isolate your problem.

*Other filters available upon request.*



Sonntek, Inc.

201-236-9300

[sonntek@aol.com](mailto:sonntek@aol.com)